

Compact Type Electrochemical Workstation ZIVE SP1

Including Internal FRA
10Volts/1Amp
Full Software Package Included

*For
Corrosion
Material Testing
Sensor/BioElectrochemistry
Battery/Fuel Cell
Super Capacitor/Solar Cell*

Electrochemical Workstation ZIVE SP1

The **ZIVE SP1** is an outstanding Potentiostat/Galvanostat/FRA offered at affordable price. This powerful model is a perfect choice for the complete DC and impedance characterization of various electrochemical applications.

The **ZIVE SP1** is equipped with a frequency response analyzer(FRA) for system as standard and it provides high performance impedance measurements over the frequency range 10uHz to 1MHz. The ZRA(zero resistance ammeter) function can measure max. 1 Amp in galvanic corrosion technique. Four(4) advanced software packages, which are categorized by application field, are provided with the system as standard. Consequently, it widens **ZIVE SP1's** flexibility.

Features

- Economical high quality Potentiostat/Galvanostat/Impedance Analyzer
- Compact size with full functions
- Smart LCD display
- Wide current ranges(1nA~1A) for various applications
- Built-in FRA : enables EIS tests by using software
- 14 EIS techniques capability including multisine
- Capable of multitude of applications
 - Corrosion, general electrochemistry, sensor, battery, fuel cell, supercapacitor, solar cell, etc.
- Current interrupt IR measurement
IR compensation(dynamic, positive feedback)
- Bipolar pulse capability
- Voltage pulse or current pulse charge/discharge test(GSM,CDMA etc.), sine wave function for ripple simulation in battery test package & pulse plating available
- High speed data sampling time
 - 50usec/sample in burst mode
 - 1msec/sample in normal mode
 - 2usec/sample in fast sweep mode
- Fast sweep mode(5000V/sec with 10mV data sampling)
- 3 measurement/control voltage ranges & 10 measurement/control current ranges
- Internal 295,000 data point storage & continuing experiment regardless of PC failure.
- Full software package included as standard
 - Corrosion test software package (CORE)
 - EIS test software package (EISe)
 - Electrochemical analysis software package (EASe)
 - Energy software package (BATE)
- Multichannel configuration available
- Free software upgrade



ZIVE SP1 front view




ZIVE SP1 rear view

Application

The **ZIVE SP1** electrochemical workstation is ideal for fundamental research in electrochemistry, development and quality assurance of new sensors, corrosion/coatings, electrode material, membrane, conducting polymer, evaluation power device research such as battery materials, fuel cells, super capacitors and solar cells.

■ General Electrochemistry




The **ZIVE SP1** is also suitable for the development of bio-research, electron transfer kinetic studies or electrochemical analysis of compounds.

■ Batteries



The system is very well adapted for researches on the cycling behavior of battery. It supports EVS(electrochemical voltage spectroscopy)/GITT/PITT test. Fast pulse capability for GSM, CDMA test is included in battery test software package. Pulse profile measurement function to check pulse shape is available. For ripple simulation test, sine wave charging/discharging is available.

■ Solar Cells



Solar cell development and production require extensive material and device testing to improve efficiency and match individual cells for panel construction. The **ZIVE SP1** is the best solution for photovoltaic cell characterization.

■ Corrosion



The system is suitable for measuring low corrosion rates and EIS test to evaluate corrosion. ZRA function is supplied for galvanic corrosion measurement.

■ Sensors



The **ZIVE SP1** can be used for sensor research using with DNA chips or screen printed electrodes. System's minimum current range is 1nA(with gain). Cyclic voltammetry, Chronoamperometry and EIS measurement can be used for this application.

■ Super Capacitors



The **ZIVE SP1** has fast potentiostat circuit with high speed data acquisition (50usec/point, burst mode). This function is well applicable to super capacitors testing. Charging/discharging capability is used for this application

Electrochemical Workstation ZIVE SP1

Fuel Cells



The ZIVE SP1 is ideal for characterizing the fuel cells and anodic/cathodic process mechanism at development and research grade. This system can be directly used for PEMFC, DMFC, and DEFC etc. Automatic current ranging potentiostatic/galvanostatic IV curve is available.

Smart Manager(SM) Software



The Smart Manager(SM) provides user defined sequential test by using sequence file, technique menu and batch file. The batch file allows the users to do a serial test by combining sequence files and/or technique files.


The SM software is easy to use and supports various electrochemical experiments including functions of system control, schedule file editor, real time graph, analysis graph, user calibration, and data file treatment etc.

Each software package's upgrade will be provided at free of charge.


Basic Techniques

Basic Technique with Standard Functions


- 1) Potentiostatic
- 2) Galvanostatic
- 3) Double step potentiostatic
- 4) Double step galvanostatic
- 5) OCP measurement
- 6) Potential sweep
- 7) Current sweep
- 8) Cyclic voltammetry
- 9) Fast potential sweep
- 10) Potentiostatic Ru measurement
- 11) Galvanostatic Ru measurement



50 usec data sampling



DC graph




Peak find module

The above functions can be used sequentially by step control function.

Sequence Editor

User can design his/her own experiment procedure by using TASK sequential routine editor.



1) Control Task Parameter

Control Mode		
constant	GSTAT	constant current control
	Crate	constant Crate control
	PSTAT	constant voltage control
	POWER	constant power control
	LOAD	constant load control
	CC-CV	constant current constant voltage control
	Crate-CV	Crate constant voltage control
	CP-CV	constant power constant voltage control
	CL-CV	constant load constant voltage control
	Id	Id control
	Is	Is control
	OCV	OCV control
Step	GSTAT	current step control
	PSTAT	potential step control
Sweep	GSTAT	current sweep control
	FAST-G	fast current sweep control
	PSTAT	potential sweep control
	FAST-P	fast potential sweep control
EIS	GSTAT	galvanostatic EIS
	PSTAT	potentiostatic EIS
	OCV	OCV EIS
	PSUEDO	pseudo galvanostatic EIS
	HFR G	galvanostatic HFR
	HFR P	potentiostatic HFR
	MisineG	galvanostatic multisine EIS
MisineP	potentiostatic multisine EIS	
Rest		rest control
ZRA		ZRA control
Loop		loop control
Pulse	Vpulse	voltage pulse control
	Ipulse	current pulse control
	GSINE	current sine wave control
	PSINE	potential sine wave control

- Constant potential, current, C-rate, power, load, OCP
- Sweep potential, current
- Fast sweep potential, current
- Staircase potential, current
- CC-CV, CP-CV, CL-CV, Crate-CV control
- Id, Is control
- EIS control
- Pulse or sinewave control
- Rest(voltage monitoring only)
- Loop(cycle) control

- 2) Cut-off(Vertex) condition
 - Time(step, test, loop, cycle)
 - Current, current density
 - Voltage
 - Capacity
 - C-rate
 - -dV
 - |dV/dt|
 - |dI/dt|
 - Aux1
 - Eoc
 - etc.

Condition-1		
Item	OP	DeltaValue
Step Time	>=	30
None		
Step End		
Step Time		
Current		
I Density		
Voltage		
Capacity		
-dV		
dt/dt		
dV/dt		
dI/dt		
Temp. (C)		
AUX1		
AUX2		
AUX3		
Test Time		
Loop Time		
Cycle Time		
Eoc		
Whr		
LCC(%)		
FCC(%)		
FCD(%)		
Power (W)		
SumQ(AHr)		
SumE(VHr)		
Loop Next		
DI Ch.		
Crate		

Cutoff Condition

- 3) Sampling Condition
 - time, |dI/dt|, |dV/dt|, |dA1/dt|, burst time
- 4) Flow View
 - This displays sequence flow at a glance.

Batch Function

User can design batch file including multiple technique files and/or sequence files. With this batch file, user can experiment several techniques/sequence in series automatically.

Index	Enable	Count	Next	Loop End	Chg	Schedule File(s)
1	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\en1.EVIS
2	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\ccv.CCV
3	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\b.LCCY
4	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\2_7u.DPE
5	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\H4.DPE
6	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\4_2V.DPE
7	<input checked="" type="checkbox"/>	1	Index-1	Next		C:\Dve Data\m\schedule\kzcv1.CCV
8	<input type="checkbox"/>	1	Next	Next		C:\Dve Data\m\schedule\kzcv.CCV


Advanced Software Packages

For a wide range of application, following software packages are provided as standard.


■ EIS Software Package (EISe)

- 1) Potentiostatic EIS
- 2) Galvanostatic EIS
- 3) Pseudo galvanostatic EIS
- 4) OCP(*1) EIS
- 5) Potentiodynamic PEIS
- 6) Galvanodynamic GEIS
- 7) Potentiodynamic HFR
- 8) Galvanodynamic HFR
- 9) Potentiostatic HFR monitor
- 10) Galvanostatic HFR monitor
- 11) Multisine potentiostatic EIS
- 12) Multisine galvanostatic EIS
- 13) Intermittent potentiostatic EIS
- 14) Intermittent galvanostatic EIS


(*1) The system measures open circuit potential before each frequency change and applies AC sine wave on this potential.




OCV EIS



Potentiodynamic PEIS



Rs, Cp & Idc vs Vdc plot




Coin Cell Intermittent PEIS
3D Nyquist Plot By ZMAN

■ Energy Software Package (BATE)

BATE software supports IR measurement.

- 1) Battery Test Technique
 - CC/CV test for cycle life test of lithium battery
 - CC/CC test for cycle life test of NiCd & NiMH battery
 - Discharging test
 - EVS (Electrochemical voltage spectroscopy)
 - Variable scan rate CV
 - Potentiostatic IV curve
 - Galvanostatic IV curve
 - Steady state CV
 - GITT
 - PITT
 - Pulse mode is available for GSM & CDMA profile.
Pulse shape profile can be measured by user's demand.



CC/CV test



EVS



Pulse shape profile monitor
(micro sec order)



Current sine wave
(charge ripple simulation)


- 2) Control Mode
 - Charge: CC, CC-CV, pulse, sine wave
 - Discharge: CC, CP, CR, pulse, sine wave

- 3) Cutoff Condition
 - Time, voltage, current, power, auxV etc.


Various battery charge/discharge test is available including pulse discharge for GSM, CDMA application.

■ Electrochemical Analysis Software Package (EASe)


- 1) Step Techniques
 - CA (Chronoamperometry)
 - CC (Chronocoulometry)
 - CP (Chronopotentiometry)
- 2) Sweep Techniques
 - LSV (Linear Sweep Voltammetry)
 - SDV (Sampled DC Voltammetry)
 - Fast CV
 - Fast LSV
- 3) Pulsed Techniques
 - DPV (Differential Pulse Voltammetry)
 - SWV (Square Wave Voltammetry)
 - DPA (Differential Pulse Amperometry)
 - NPV (Normal Pulsed Voltammetry)
 - RNPV (Reverse Normal Pulse Voltammetry)
 - DNPV (Differential Normal Pulse Voltammetry)



50 usec sampling



sdv_charging/sdv_discharging




NPV & RNPV overlay


■ Corrosion Software Package (CORE)

Corrosion technique supports IR compensation.


- 1) Tafel (Tafel Experiment)
- 2) Rp (Polarization Resistance)
- 3) RpEc Trend
- 4) PDYN (Potentiodynamic)
- 5) CYPOL (Cyclic Polarization Resistance)
- 6) GDYN (Galvanodynamic)
- 7) Reactivation
- 8) Galvanic Corrosion
- 9) Potentiostatic ECN
- 10) Galvanostatic ECN
- 11) ZRA mode ECN



Tafel plot



Polarization resistance



Cyclic polarization resistance

Control & Real Time Graph


Smart Manager Program provides virtual control panel for control & data acquisition with real time graph.



Multichannel real time graph

User can control and monitor in this control panel and he/she can monitor data in VOI (Value Of Interest) window and channel status in one window. Real time graph's X Y axis format will be changed per technique automatically. It can be defined by user's demand per technique. For experiment using sequence file or batch file, user can designate X,Y axis parameter on three different real time graph.

Each real time graph format can be also selected. Real time graph and VOI will be changed depending on DC test or impedance test automatically. Virtual control panel always displays the graph for recent test result. For impedance measurement, wave monitor will be displayed on real time graph to check wave's quality. This monitor can be switched to Lissajous (I vs. E) plot.




Single channel control panel



Multichannel real time graph

Strip Chart







Strip chart recorder function allows you to monitor DC curve in real time. You can monitor 2 Y axis data such as voltage, current, auxV, power, capacity etc. in real time


Graph



Smart Manager's graph function is to simplify the operation. There are 3 kinds of graph per each experiment. You can change X, Y1, Y2, Y3, Y4 axis parameter as you want. Each graph provides shortcut buttons. When you click these buttons, the format of the graph will be changed accordingly.




In DC and Cycle graph, whenever you click  or , the parameters which are related to current such as current, capacity, energy, power, load, etc., are changed into calculated specific value or density value, respectively.

 : value divided by weight

 : value divided by active area




1) DC Graph

- For general data display
- 4 shortcut buttons: I vs. V, E vs. LogI, V, I vs. time, V vs. Q
- Graph parameters: Time, Eref, I, Eoc, Id, Aux1, LogI, Load, ChQ, DchQ, ChQs, DchQs, Ch P, Dch P, Ch-Wh, Dch-Wh, Sum Wh, Sum Q, Sum |Q|, |Q|, Rp, dQ/dV



2) EIS Graph

- For EIS data display
- 3 shortcut buttons: Nyquist plot, Bode plot, Cs vs. frequency
- Graph parameters: Frequency, Zre, -Zim, Zmag, Zph, Y, Yimg, Y, |Y|, Yph, LogZ, LogY, Rs(R-C), Cs(R-C), Rp(R|C), Cp(R|C), Rs(R-L), Ls(R-L), Q(R-L), time, Vdc, Idc, Aux1





3D Bode Plot by ZMAN


Technique used: Potentiodynamic Impedance Measurement
By Using A Corrosion Cell

3) BAT Graph


- For battery cycle data display
- 3 shortcut buttons: cycle capacity, cycle average, Log(cycle No) vs. depth of discharge plot.
- Graph parameters: cycle number, Ch Q, Dch Q, Sum Q, Coulomb Eff, Ch-Wh, Dch-Wh, Sum Wh, Energy Eff, MinV, MaxV, ChQs, DchQ, ChVavg, DchVavg, Vavg




■ Data Export to ASCII & Excel File



Selectable between 'Convert data on graph only' and 'Convert selected file(s)'



3D Bode plot for series measurement




Parameter plot

Data Analysis Software


ZIVE data file can be used for analysis by using external IVMAN™ software for DC analysis, IVMAN DA™ software for battery data analysis, IVMAN PA™ software for photo-voltaic cell data analysis and ZMAN™ software for EIS data analysis without license.

ZMAN™ EIS Data Analysis Software


- Model simulation and fitting
- 2D- and 3D-Bode- and Nyquist plots
- Automatic equivalent circuit model search function
- Project concept to handle multiple EIS data analysis
- Parameter plot from fitted elements value
- Compatible with data format from Zahner, Gamry, Ametek etc. (License code is needed.)
- Various weighting algorithm
- Model library and user model
- KK plot
- Batch fitting for project data
- Impedance parameter simulation
- Interpolate bad data
- Black-Nichols plot
- 3D graph setting option
- Improved model editor
- Application model library for automatic searching
- Parameter simulation of model
- Genetic algorithm option for initial guessing
- Automatic initial guessing
- Trace movie function on fitting
- Free for ZIVE's data format(*.seo, *.wis) analysis (No license code required.)
- Circle fitting
- Data editing available (insert, delete, edit)
- Add/subtract element parameters
- Add/subtract model parameters
- Impedance, Z in polar, admittance, Y in Polar, modulus, M in polar, dielectric constant, E in polar. data display
- Empty cell capacitance calculation
- Find file function
- Data replacement by formula function
- Cursor data display
- Model finding result automatic sorting by Chi square value
- R, C R, L R, Q preview & graphic
- ZHIT function
- Mott-Schottky analysis
- Donor density vs. Vfb graph
- C vs. voltage graph




Importing 3rd parties ASCII data file




Project manager with data preview





Cursor data display





Data replacement by formula function




Fitting display



Model editor & model library



Automatic model searching



LEVM fitting

Electrochemical Workstation ZIVE SP1

IVMAN™ DC Data Analysis Software




- IVMAN™ software package consists of
 - IVMAN software
 - IVMAN utilities
 - IVMAN differential analysis software
 - IVMAN photo voltaic cell analysis.
 - IVMAN Tafel analysis
 - IVMAN extractor
 - IVMAN peak find module




IVMAN DA™ Battery Test Data Analysis Software


- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. Capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph(dV/dQ vs. Q)




Parameter simulation




Element add/subtraction




Finding data file menu




Circular fitting




2D Nyquist plot




2D Bode plot




Rp,Cp vs frequency (R/C)




Empty cell capacitance




Donor density vs. Vfb graph and analysis




KK consistency



Mott-Schottky analysis window




C/R-V graph




Measured data

V vs. Q



Cycle graph

dQ/dV vs. V




dV/dQ vs. Q

Export ASCII file



IVMAN™ Photovoltaic Cell Analysis




- Automatic analysis of parameters
 - open circuit voltage, open circuit current, max. power, efficiency
 - photo induced current, diode quality factor, series resistance, etc.




IVMAN TA™ Tafel Analysis

- Simple Tafel calculation



IVMAN EX™ Extractor

- Extracting data by cycle number or step
- Exporting ASCII file








IVMAN™ Main Software

- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis
- Polarization resistance fitting
- 3D graph
- Find peak function
- Interpolation, differentiation, integration etc.
- Reporting function

IVMAN PF™ Peak Find Module

- Independent peak finding software



Optional Accessories

- Power Booster
 - for high voltage/high current application
 - modular type design
 - EIS capability
 - sine wave simulation available




- Battery Jig & Coin Cell Jig
 - for cylindrical cell or coin cell
 - 4 probe type



Electrochemical Workstation ZIVE SP1

- Pouch Cell Jig
 - contact type
 - a) pull-down contact type with adjustable contact probe's width
 - b) banana connector for cell cable connection
 - 4 contact point type(Kelvin probe)



- Through-Plane Conductivity Test Jig
 - for through plane conductivity measurement
 - 2 probe type




- Coin Cell Holder
 - D-SUB connector type



- Single Cell Hardware Fixture
 - for PEMFC and DMFC
 - max. temp. : 120°C or 180°C
 - active area : 5, 9, 25, 50, 100cm²
 - MEA is not included.



- Membrane Conductivity Cell
 - for 5, 9 and 25cm² fuel cell hardware fixture
 - material : PEEK(cell body), platinum(wire)
 - operating temp. : up to 130°C



- Battery Test Cell
 - test cell for 2 or 3 electrode measurement
 - test cell for aqueous electrolyte
 - test cell for pressure monitoring
 - test cell for gas diffusion electrode
 - test cell for time-resolved gas analysis
 - test cell for optical and X-ray characterization in the reflective mode
 - test cell for high-throughput testing, etc.

- Universal Electrode Holder
 - electrode and glass vial are not included.

- Faraday Cage
 - size : 300 x 300 x 398mm (WxDxH)



ECC-Ref



ECC-Opto



ECC-Press



PAT-Single Stand

- Cell Kit



Corrosion Cell Kit



Flat Cell Kit

- Electrochemical Dilatometer
 - for the measurement of height changes in aprotic as well as aqueous electrolytes
 - ECD-3 : displacement resolution ≤ 50 nm
 - ECD-3-nano : displacement resolution ≤ 5 nm



Plate Test Cell



Plate Test Cell



H-Type Cell



H-Type Cell

Specification

Main System	
PC communication	USB2.0 high speed
Line voltage	100~240VAC, 50/60Hz, 1Amp
Power adapter	24V 2.5Amp
Size/weight	160x330x81mm(WxDxH) / 2.05Kg
Max. output power	15Watt

System	
Cell cable	1 meter shielded type(standard) working, reference, counter, working sense, Auxiliary V
Control DAC	DSP with FPGA 2x16bit DAC(50MHz) for bias & scan
Data acquisition ADC	2x16bit ADCs(500kHz) for voltage, current 1x16bit ADCs(250kHz) for auxiliary reading
Calibration	Automatic
Filter selection	4ea(5Hz, 1kHz, 500kHz, 5MHz)
Scan rate	0~200V/sec in common mode 0~5000V/sec in fast mode
Max. channel No.	8 channels via USB connection
Internal data memory	295,000 points

Power Amplifier(CE)	
Power	12Watt (12V@1A)
Compliance voltage	±12V
Max. current	±1A
Control speed selection	4ea
Bandwidth	2MHz
Slew rate	10V/usec

Potentiostat Mode (voltage control)	
Voltage control	
Control voltage range	±10V, ±1V, ±100mV
Voltage resolution	16 bit per each range
Voltage accuracy	±1mV ±0.05% of setting(gain x1)
Max. scan range	±10V vs. ref. E
Current measurement	
Current range	10 ranges(auto/manual setting) 100nA ~ 1A 1nA & 10nA with gain
Current resolution	16 bit 30uA, 3uA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA (300fA, 30fA with gain)
Current accuracy	±10pA ±0.1% f.s.(gain x1)>100nA

Galvanostat Mode (current control)	
Current control	
Control current range	max. ±1A ± full scale depending on selected range
Current resolution	16 bit 30uA, 3uA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA (300fA, 30fA with gain)
Current accuracy	±10pA ±0.1% f.s.(gain x1)>100nA f.s.
Voltage measurement	
Voltage range	10V, 1V, 100mV
Voltage resolution	16 bit 0.3mV, 30uV, 3uV
Voltage accuracy	±1mV ±0.05% of reading(gain x1)

Electrometer	
Max. input voltage	±10V
Input impedance	2x10 ¹³ Ω 4.5pF
Bandwidth	>22MHz
CMRR	>114dB

EIS(Internal FRA) for System	
Frequency range	10uHz~1MHz
Frequency accuracy	0.01%
Frequency resolution	5000/decade
Amplitude	0.1mV~5V rms (Potentiostatic) 0.1~70% f.s. (Galvanostatic)
Mode	Static EIS: Potentiostatic, Galvanostatic, Pseudogalvanostatic, OCP Dynamic EIS: Potentiodynamic, Galvanodynamic Fixed frequency impedance: Potentiostatic, Galvanostatic, Potentiodynamic, Galvanodynamic Multisine EIS: Potentiostatic, Galvanostatic Intermittent PEIS/GEIS

Interfaces for System	
Auxiliary port	
Auxiliary voltage input	1 analog input: ±10V
Zero resistance ammeter	100nA~1A ranges
External booster interface	Via booster I/F cable
Sync terminal	For channel synchronizing

Smart LCD Display	
DC mode	Control value, E value, I value E range, I range
EIS mode	Frequency, Magnitude, Phase E range, I range
Operation status	Mode: PST, GST, ZRA, EIS, CC, CV, CP, CR Status: Cell On, Run, Error

Software	
Max. step per experiment	1000
Shutdown safety limits	Voltage, current, power, AuxV etc.
Max. sampling rate	20kHz(50usec) in burst mode 500kHz(2usec) in fast sweep mode
Min. sampling time	Unlimited
Sampling condition	Time, dv/dt, di/dt, etc.

PC Requirement	
Operating system	WindowsXP SP3/7/8/10(32bit/64bit OS)
PC specification	Pentium4, RAM 1GB or higher
Display	1600x900 high color or higher
USB	High speed 2.0

General	
Dummy cell	One external dummy cell included
Impedance analysis S/W	ZMAN™ software
DC data analysis S/W	IVMAN™ software package

The specifications are subject to change without notice.
Windows is a registered trade mark of Microsoft Corporation.

Designed by

ZIVE LAB
www.zivelab.com



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

eDAQ Pty Ltd
6 Doig Avenue
Denistone East
NSW 2112 Australia
+61 2 98078855
e-mail info@edaq.com
web www.edaq.com



ISO 9000 & ISO 14000 Qualified



CE Certified