



C⁴D Headstage for Capillary Electrophoresis (ET120)



C⁴D Headstage for CE

- Compatible with most capillary electrophoresis instruments
- Use with silica or plastic capillary (350 to 385 µm OD)
- Compatible with all eDAQ C⁴D units
- Internal electrodes never require maintenance or cleaning

Description

The ET120 headstage can be connected to any eDAQ C⁴D unit to make a contactless conductivity detector (C⁴D) for capillary electrophoresis. The capillary tube (350 to 385 µm outer diameter) is pushed through this headstage, where a pair of internal ring electrodes are positioned. A high frequency AC signal is passed between the electrodes, and the conductivity of the solution in the capillary affects the received AC amplitude. The C⁴D unit generates the excitation waveform, and provides signal filtering and offset.

Installation is easy: simply slide the capillary through the headstage. There is no need for windowing (scratching off the capillary's polyimide coating) as with optical detectors.

There is a model of the ET120 available for capillaries with 150 µm outer diameter.

Compatibility

Can be used with any eDAQ C⁴D unit.

Specifications

Mechanical	
Guide Tube:	400µm ID, 630µm OD
Suitable Capillary:	350 to 385 µm OD
Electrode length:	2mm
Electrode separation:	2mm
Size:	25mm x 34mm x 13mm
Weight:	20grams
Electrical	
Output voltage:	-0.100 to +5.00Volts
Linear Range:	+0.100 to 2.500Volts

The headstage adapter cable is included with the C⁴D unit, not with the ET120. The correct adapter depends on which C⁴D unit is being used:

EC1210 headstage adapter is used with:

- ER815 C⁴D Detector.
- ER825 C⁴D Multi channel detector.

EC1208 headstage adapter is used with:

- ER225 C⁴D Data System.
- EA120 C⁴D Amplifier.
- ER125 C⁴D Detector.

The ET120 is compatible with most CE systems, including Agilent, Beckman Coulter and PrinCE instruments.

Applications

Contactless conductivity detection can be used for virtually all charged species: inorganic anions and cations, as well as organic ions, such as carboxylic acids, amines, amino acids, peptides, proteins, DNA fragments, antibiotics and many other pharmaceutical compounds.

Excitation Level:	1 to 100 Volts peak to peak
Excitation Frequency:	50kHz to 1200kHz +/-1db
Gain x1 =	100mV/µA
Gain x5 =	500mV/µA
Frequency response:	0-10Hz (3dB point)
Noise:	< 3 µV RMS 0-10Hz measured over 1 sec
Drift:	< 3µV/°C measured with no capillary
Power supply:	+/-5Volts +/-5%
<i>eDAQ Pty Ltd reserves the right to alter these specifications at any time.</i>	

www.eDAQ.com

E-mail: info@edaq.com

e-corder is a registered trademark of eDAQ Pty Ltd.
All other trademarks are the property of their respective owners.

Document Number: M-ET120-1405

Copyright © eDAQ 2014