



ET014 EChem Electrode Kit

Description

The ET014 EChem Electrode Kit includes:

- three working electrodes
 - platinum
 - gold
 - glassy carbonin PEEK (polyetheretherketone) body, with dimensions of 3 mm OD x 67 mm length, and 1 mm diameter active area
- leakless Ag/AgCl reference electrode
- platinum coated auxiliary electrode
- two 3 mL reaction vials with stirrer vanes
- cell stand
- Helping Hands cable holder



Ordering Information

ET014 EChem Electrode Kit

For use with:

EA160 Potentiostat
EA161 Potentiostat
EA162 Picostat
EA164 Quadstat

www.eDAQ.com

E-mail: info@edaq.com

e-corder is a registered trademark of eDAQ Pty Ltd.
All other trademarks are the property of their respective owners. TT7/03

Document Number: MET014-0409
Copyright © eDAQ 2009