

ET014 EChem Electrode Kit

An electrode kit for conducting electrochemical experiments. It is suitable for use with both the EA160 Potentiostat and EA162 Picostat.

The Kit includes:

- three different working electrodes encased in PEEK (3 mm OD x 67 mm)
 - platinum (1 mm diameter active area)
 - gold (1 mm diameter active area)
 - glassy carbon (1.5 mm diameter active area)
- a 67 mm Ag/AgCl reference electrode
- a 67 mm platinum auxiliary electrode
- two 3 mL reaction vials
- a cell stand



ET014 EChem Electrode Kit