



ET1121 Membrane Housing Kit

Description

The ET1121 Membrane Housing Kit includes:

- 6 replacement Teflon membranes
- a bottle of electrolyte solution
- a bulb pipette
- 5 polishing pads



Ordering Information

ET1121 Membrane Housing Kit (for ET1120)

For use with:

ET1120 Micro-Oxygen Electrode

www.eDAQ.com

E-mail: info@edaq.com

e-corder is a registered trademark of eDAQ Pty Ltd.
All other trademarks are the property of their respective owners. 117/03

Document Number: UET1121-0703
Copyright © eDAQ 2003